


Assessing In Vitro Drug Release from Multivesicular Liposome: Comparison of Reverse Dialysis and Rotary Shaking Methods

Soumyarwit Manna^{1,2}, Peter Petrochenko¹, Yong Wu², Bonhye Koo^{1,2}, Ke Ren¹, Lynn Chen^{1,2}, Yixuan Dong¹, Xiaoming Xu¹, Stephanie Choi¹, Darby Kozak¹, Yan Wang¹, Jiwen Zheng²

¹ Division of Therapeutic Performance, Office of Research Standards, Office of Generic Drugs, Center for Drug Evaluation and Research, ² Division of Biology, Chemistry and Materials Science, Office of Science and Engineering Laboratories, Center for Devices and Radiological Health, Food and Drug Administration, Silver Spring, MD 20993

Email : soumyarwit.manna@fda.hhs.gov



PURPOSE

- Bupivacaine (BPV), an amide type local anesthetic, has a relatively short half-life ($t_{1/2-BPV}$) ~ 2.7 hours. BPV, when encapsulated in a multivesicular liposome (MVL) formulation, exhibits sustained release characteristics
- Multivesicular liposomes (MVLs) consist of a non-lamellar honeycomb structure with non-concentric aqueous chambers
- Significantly increased elimination $t_{1/2-BPV}$ ~ 34 hours is observed, which is likely influenced by various physicochemical parameters of the MVL such as:
 - MVL vesicle morphology and size (both external and internal physical dimension)
 - Rate of drug partitioning between hydrophilic (aqueous phase) and lipophilic (lipid bilayers) components
- It is *hypothesized* that alteration of MVL morphology over time influences the in vitro drug release rate

OBJECTIVES

- Characterize the vesicle size and morphology of the BPV-MVL complex formulation
- Conduct in vitro drug release in 50 mM PBS with and without Human Serum Albumin (HSA) over a period of 7 days - expected complete drug release profile (5 times the $t_{1/2-BPV}$)
- Analyze the effect of different parameters such as a) agitation, b) temperature, c) dissolution medium, and d) dissolution set-up on the in vitro drug release rate

METHODS

- MVL vesicle size – Laser diffraction (Malvern Mastersizer 3000)
- MVL vesicle morphology – Cryo-scanning electron microscopy (Tescan MIRA 3)
- In vitro dissolution set-up and drug release detection techniques:
- Rotary shaker mechanism (Thermo Scientific™ Tube Revolver / Rotator)

Dissolution media	50 mM PBS (pH 7) with varied concentrations of HSA		
Temperature	37°C	Rotation	12 rpm
		Dilution of model drug	17 x

HPLC Column	Agilent ZORBAX SB-CN – 4.6x150 mm, 5µm		
Mobile phase	40% Acetonitrile + 0.01% Trifluoroacetic acid; pH 2.8		
Flow rate	1 mL/min	Temperature	30°C
λ_{BPV}	263 nm	Retention time	~ 2.1 min

- Reverse dialysis (USP II apparatus – Teledyne Hanson Research)

Dissolution media	50 mM PBS (pH 7)		
Temperature	37°C and 40°C	Rotation	120 rpm and 240 rpm

Dialysis cartridges	10-100 kDa (MWCO) conditioned with SDS
Drug detection	In situ UV-Vis fiber optic probe
Dilution of model drug	17x


Rotary shaker device




Reverse Dialysis - USP II apparatus




RESULTS: MVL Vesicle Size and Morphology



RESULTS: In vitro Drug Release – Rotary Shaker - HPLC



RESULTS: In vitro Drug Release – Reverse Dialysis USP II-in situ UV



Normal dialysis of BPV standards using 100 kDa MWCO cartridges shows a) better detection of initial release of BPV (~2 h) and b) 100% release of BPV

Acknowledgments: The authors would like to thank the FDA Advanced Characterization Facility (ACF) and CDRL/Office of Science and Engineering Laboratories/Division of Biology, Chemistry and Material Science for instrument use.
Dr. Manna was supported in part by an appointment to the Oak Ridge Institute for Science and Education (ORISE) Research Participation Program at the Center for Drug Evaluation and Research administered by the ORISE through an agreement between the U. S. Department of Energy and FDA/CDER.

Disclaimer: The views expressed in this poster do not necessarily reflect the official policies of the Department of Health and Human Services; nor does any mention of trade names, commercial practices, or organization imply endorsement by the United States Government.

CONCLUSIONS

- BPV MVLs are monodisperse spherical particles with complex internal compartment morphology
 - In vitro release – **Rotary shaker** → Orientation of agitation and presence of HSA influence BPV release
 - In vitro release – **Reverse dialysis – USP II coupled with in-situ UV fiber optic** → Discern both an initial diffusion burst release (~ 2 – 4 hr) followed by secondary release likely triggered by physical degradation of MVLs. It also provides the advantage of continuous monitoring of the BPV release profile.
- Temperature ↑ initial release ↑ -- Agitation speed ↑ Total BPV release duration ↓