

Comparing Nasal Suspension Products Using Realistic *In Vitro* Test Methods

Michael Hindle and Mandana Azimi
Department of Pharmaceutics
Virginia Commonwealth University


Goals

Development and evaluation of methods to characterize nasal spray products using realistic nasal airway models as more clinically relevant *in vitro* tools:

- Patient use variability
- Inter-subject variability
- Product variability

Nasal drug delivery

- ❑ Can be used for local or systemic delivery
- ❑ Metered dose nasal sprays are the most commonly used devices
 - Relatively poor delivery efficiency to the site of action in the middle passages
- ❑ Drug delivery efficiency depends on:
 - Patient use
 - Nasal geometry
 - Formulation and device combination




Nasal spray product characterization

- ❑ Nasal drug delivery efficiency and assessments of bioequivalence currently use *in vitro* characterization methods that focus on the spray plume and droplet size characteristics of the nasal spray.
- ❑ Statistical differences spray plume properties may not lead to changes in nasal drug deposition¹ which questions the clinical value of the current *in vitro* tests.
- ❑ Assessments of nasal drug delivery efficiency and bioequivalence may be aided by the use of more clinically relevant *in vitro* testing using **physically realistic nasal airway models** combined with **simulated patient use parameters**.

¹Suman et. al, *J Aerosol Med*, 2006

Clinically relevant *in vitro* nasal testing



Regional drug deposition measured on:

i) Nasal spray device


iii) **Middle passages + nasopharynx**

ii) Anterior nose region + drip

iv) Throat + filter

Innovator product: Nasonex[®] Nasal Spray (mometasone furoate monohydrate) 50 µg/100 µl

Nasal geometry: VCU nasal model 1




Data set

Guilmette data, MRI scan of an individual - VCU Model 1

Dh, nostril and nasopharynx	12.1 mm, 5.9 mm
Surface area (SA)	8024.2 mm ²
Volume (V)	10832mm ³
SA/V	0.7 mm ⁻¹
SA of the nasal valve	1156 mm ²
Anterior nose volume	3.2 ml

Patient-use variables and DoE design


Head angle: 30° or 50°




Position: 9 or 5 mm




Actuation force: 4.5 or 7.5 kg



Timing: D or E




Middle passage deposition: Mometasone innovator product




- Nasal deposition varied significantly with changing patient use factors
- Significant main effect of nasal spray position within the nostril
- Significant interaction between inhalation timing & head angle

Inter-subject variability




Data set	Guilmette data, MRI scan of an individual - VCU Model 1	VCU Medical Center, CT scan of an individual - VCU Model 2
Dh, nostril, nasopharynx	12.1 mm, 5.9 mm	10.6 mm, 4.5 mm
Surface area (SA)	8024.2 mm ²	6802.3 mm ²
Volume (V)	10832mm ³	5118 mm ³
SA/V	0.7 mm ⁻¹	1.3 mm ⁻¹
SA of the nasal valve	1156 mm ²	1493 mm ²
Anterior nose volume	3.2 ml	2.2 ml

Middle passage deposition: Mometasone innovator product



- Lower impact of patient use factors on nasal deposition in VCU Model 2
- Significant effect of inhalation timing & head angle

Middle passage deposition: Mometasone innovator product




○ Higher middle passage deposition in VCU model 2 compared to VCU model 1


Mean middle passage drug deposition (% recovered dose) and standard deviation, n= 4. *p<0.05 (student t-test).

Comparison of *in vitro* and *in vivo* deposition data for Nasonex[®] nasal spray product

Study 1 (Leach et. al , 2015)




Study 2 (Shah et. al, 2014)





- Mean (SD) middle passage deposition in VCU models 1 and 2 were 34.0 (13.6)% and 64.8 (12.4)%, respectively across all test conditions.

Regional nasal drug deposition of innovator & generic nasal spray products

- ❑ Formulation and device
 - Mometasone furoate: innovator vs “in house” (University of Bath) nasal spray
 - Fluticasone: innovator vs generic nasal spray
- ❑ Nasal geometry: VCU models 1 & 2
- ❑ Patient-use conditions: low, intermediate and high middle passage deposition



Middle passage deposition: Mometasone innovator and “in house” products




- **No statistical difference in the middle passage drug deposition for the two nasal spray products at each respective level**

Mean regional deposition (% recovered dose) and standard deviation, n= 4.

Middle passage deposition:

Fluticasone innovator and generic products




- **No statistical difference in the middle passage drug deposition for the two nasal spray products at each respective level**

Mean regional deposition (% recovered dose) and standard deviation, n= 4.

Regional nasal drug deposition of mometasone furoate nasal sprays with varying spray plume characteristics

- ❑ Formulations:
 - “In house” mometasone furoate nasal spray: Batch A
 - “In house” mometasone furoate nasal spray: Batch B
 - Mometasone innovator nasal spray: expiry Feb 2007
 - Mometasone innovator nasal spray: expiry Oct 2015 (in-date)
- ❑ Nasal Geometry: VCU model 1
- ❑ Patient-use conditions: Intermediate level
- ❑ FDA recommended *in vitro* tests for establishing BE for nasal spray products
 - Droplet size distribution (median volume diameter, D_{v50}) measured by Malvern Spraytec
 - Spray pattern (smallest diameter on spray plume image, D_{min}) measured by SprayVIEW™



Mometasone furoate nasal spray middle passage deposition

Device	D _{v50} (μm)	D _{min} (mm)	MP (%)
Batch A	35.6 (0.5)	28.3 (0.6) ^a	31.0 (8.5) ^a
Innovator 2007	36.3 (0.8)	23.8 (0.5) ^{a, b}	38.8 (3.8) ^a
Batch B	38.7 (3.0)	21.3 (0.2)	49.6 (8.8)
Innovator 2015	47.0 (1.0)	20.5 (0.2)	44.6 (7.6)

Mean (SD) for D_{v50} and D_{min} for mometasone furoate nasal spray products and middle passage drug deposition (MP), expressed as % recovered dose, n = 3 and 4.

P<0.05, ^a compared to Batch B, ^b compared to Batch A

- **Lower middle passage drug delivery observed with larger plume diameters (D_{min})**

Mometasone furoate nasal spray middle passage deposition

Device	D _{v50} (μm)	D _{min} (mm)	MP (%)
Batch A	35.6 (0.5)	28.3 (0.6)	31.0 (8.5)
Innovator 2007	36.3 (0.8)	23.8 (0.5)	38.8 (3.8)
Batch B	38.7 (3.0) ^c	21.3 (0.2)	49.6 (8.8)
Innovator 2015	47.0 (1.0)	20.5 (0.2)	44.6 (7.6)

Mean (SD) for D_{v50} and D_{min} for mometasone furoate nasal spray products and middle passage drug deposition (MP), expressed as % recovered dose, n = 3 and 4.

P<0.05, ^c compared to Innovator 2015

- Differences in plume geometries appeared to be more critical than droplet size in determining the deposition fate of the nasal spray

Conclusions and Future Studies

- ❑ A realistic *in vitro* test was developed for nasal sprays incorporating realistic patient use conditions and airway nasal models
- ❑ Patient use factors and nasal model geometry were observed to have significant effects on middle passage drug delivery
- ❑ Innovator and generic nasal spray products were observed to have similar regional *in vitro* nasal deposition profiles
- Studies required to further understand relationship between *in vitro* spray characterization test metrics and nasal spray deposition.
- Develop nasal geometries that capture mean and high/low range of nasal spray deposition across the adult population.

Acknowledgements

Dr. Worth Longest, Virginia Commonwealth University
Drs. Guenther Hochhaus and Jurgen Bulitta, University of Florida,
Drs. Robert Price and Jag Shur, University of Bath
Drs. Renish Delvadia, Bhawana Saluja, Mohammad Absar, FDA.

Partial funding was provided by Contract # HHSF223201310220C,
from the Department of Health and Human Services, Food and Drug
Administration.

Views expressed in this presentation do not necessarily reflect the
official policies of the Department of Health and Human Services; nor
does any mention of trade names, commercial practices or
organizations imply endorsement by the United States Government