

SEP 26, 2022



TEXAS

The University of Texas at Austin

MELT-EXTRUDED DEXAMETHASONE OPHTHALMIC IMPLANTS: PROCESS, STRUCTURE AND IN VITRO DRUG RELEASE


This work was supported by the Broad Agency Announcement (BAA) Contract # 75F40120C00198 from the U.S. Food and Drug Administration (FDA). The content reflects the views of the authors and should not be construed to present FDA's views or policies.

FENG ZHANG

The University of Texas at Austin, College of Pharmacy
Department of Molecular Pharmaceutics & Drug Delivery

12th American Drug Delivery & Formulation Summit
San Diego, CA


Dexamethasone (DEX) Ophthalmic Implant



| Component | % wt | wt. (mcg)/implant |
|---------------------------------------|------|-------------------|
| Dexamethasone, Form B (micronized) | 60 | 700 |
| 50:50 PLGA ester with acid end group | 30 | 350 |
| 50:50 PLGA ester with ester end group | 10 | 116 |


Dexamethasone

- Anti-inflammatory corticosteroid
- Poor aqueous solubility: 90 $\mu\text{g}/\text{mL}$
- LogP: 1.83
- Non-ionizable
- $T_m = 265\text{ }^\circ\text{C}$
- Two crystalline forms: A and B





DEX degrades upon melting

- TGA shows significant degradation of DEX occurs once it melts
- No weight loss below 150°C
- Melt temp for melt extrusion of ophthalmic implant: 100-110°C



Little interaction between DEX and PLGA

DEX does not depress PLGA T_g




Gasmi H, Siepmann F, Hamoudi MC, Danede F, Verin J, Willart JF, et al. Towards a better understanding of the different release phases from PLGA microparticles: Dexamethasone-loaded systems. *Int J Pharm.* 2016;514(1):189-99.


PLGA thermally stable below 140°C

- PLGA can degrade by hydrolysis and thermal mechanisms
- TGA shows thermal degradation is not significant below 200°C
- RG502 and RG502 H PLGAs were demonstrated to be stable below 140°C during HME
- Adsorbed moisture accelerates degradation of RG502 H > 30% RH

MW change after melt extrusion




MW change after heat treatment (140°C, 30 minutes)




In-house DEX implant manufacturing process

Schematic




In-house ocular implant manufacturing process

Video



Accurate diameter control achieved with 2-stage feeding


Steady State Diameter
(n = 1474)




Structural analysis of DEX implant

SEM and MicroCT reveal irregular surface and 6% internal porosity

SEM profile




MicroCT profile



**MicroCT
Cross-section**

Internal voids
seen as dark spots

2-step extrusion process for implant structure control



2-stage feeding process to achieve consistent feed rate (8-10 g/hr)

Possible solutions to improve feeding performance of Haake force feeder

1. Volumetric/gravimetric feeder → granules bridge over screws, need stuffer feeder
2. Force feeder with more powerful motor
3. **Combination of volumetric (Brabender) and stuffer feeding (Haake)**


Stage 1: Brabender twin-screw volumetric feeder

Stage 2: Haake force feeder (water-cooled jacket to just above dew point)



Puller for controlled drawdown to target diameter


Purpose-built unit for fine adjustment of drawdown parameters



Long residence time observed due to need for slow and steady feeding (8-10 g/hr)



- Mean residence time = 33 min
 - Used tracer dye and laboratory camera to detect color change
 - In-house software developed for data acquisition and analysis
- Assay confirms DEX stability and uniformity
 - Potency: $100.1 \pm 0.2\%$ SD (N=6)
- GPC data shows 7% reduction in M_n

Residence Time Distribution Analysis
(Second Extrusion, N=94440)



Steady-state torque correlates with implant porosity

Steady-state torque can be thought of as extruder back-pressure



Increasing steady-state torque

DEX Ophthalmic Implant – Process and Structure


- Low glass transition temperature of PLGA enables the extrusion at relatively low temperature to maintain the chemical stabilities of dexamethasone
- Accurate diameter control was achieved with a two-step extrusion process using two-stage feeding design and a purpose-built downstream puller
- Various features of implant structure identified are:
 - diameter
 - pore size/distribution
 - surface roughness
- Preliminary data indicates that:
 - DEX does not dissolve in PLGA melt during extrusion and is dispersed in PLGA matrix as crystalline particles
 - DEX is stable during the current melt extrusion process

Dissolution testing of dexamethasone implant

Excerpt from patent

FIGS. 2 and 10 also demonstrate that after 28 days in vivo in rabbit eyes, or in vitro in a saline solution at 37°C., respectively, almost all of the active agent has been released from the implants. Furthermore, FIGS. 2 and 10 **show that the active agent release profiles for the extruded implants in vivo (from the time of implantation) and in vitro (from the time of placement into a saline solution at 37° C.) are substantially similar and follow approximately a sigmoidal curve, releasing substantially all of the active agent over 28 days.**

in vitro / in vivo release from Ozurdex patent data




Replotted based on the data (Figs 2 and 10) in the patent

In vitro release testing in normal saline of DEX implants

Tri-phasic release profile aligns with published data


Dexamethasone Intravitreal Implant Dissolution

37°C, 30 mL, Normal Saline, N=6 (mean ± S.D.)



Moisture Pickup during Dissolution

37°C, 30 mL, Normal Saline, N=6 (mean ± S.D.)




In vitro release testing in normal saline of DEX implants

Limited drug release in first week despite substantial changes to the implant


GPC - Change in PLGA Molecular Weight

37°C, 30 mL, Normal Saline, N=3 (mean ± S.D.)




Change in Media pH

37°C, 30 mL, Normal Saline, N=6 (mean ± S.D.)











In vitro release testing in normal saline of DEX implants

SEM shows significant structural changes to implant after day 7



In vitro release testing in normal saline of DEX implants


SEM shows significant structural changes to implant after day 7



Similar tri-phasic release profile with slower release rate in phosphate buffer pH 7.4

Dexamethasone Intravitreal Implant Dissolution

37°C, 100 RPM, 30 mL, N=3



DEX Ophthalmic Implant – Drug Release

- Tri-phasic profile derives from both implant structure and physicochemical interactions between DEX and PLGA
- Limited solubility of DEX in PLGA results in DEX crystals uniformly dispersed throughout the PLGA matrix after melt extrusion
- Burst release at the onset was not observed, due to inaccessibility of DEX crystals coated with PLGA on implant surface
- In the initial lag phase, significant PLGA undergoes significant hydrolysis even though less than 10% drug is released
- Dexamethasone release is controlled by the erosion of the implant matrix

Acknowledgements



- Mark Costello, Graduate Student
- Dr. Beibei Chen, Postdoc

U.S. Food & Drug Administration,
Office of Research and Standards,
Office of Generic Drugs, CDER

- Dr. Bin Qin, Staff Fellow
- Dr. Yan Wang, Team Lead